# **VISUAL ACUITY SERIES**



#### **DESCRIPTION**

## **Distance Eye Charts**

This Visual Acuity Chart Collection was designed by an optometrist, Anthony Hanks OD.

#### **Testing Distance:**

The standard testing distance for visual acuity (VA) is 6 metres (or 20 feet), however sometimes there are space limitations in situations like a general practitioner's medical practice, or during a home visit to a patients with mobility issues. On these occasions a chart for a shorter testing distance of 3 metres (10 feet) is useful.

#### Note re the "Bottom Line":


Most charts have "normal vision" (6/6 or 20/20) as the bottom line, since most patients like to be able to read to the last line. However, some verisons here also have one line better than normal (6/5) for practioners who prefer this.

#### **High Resolution**:

All charts in this series have been produced at a high resolution of 400 dpi.

#### **Reference Bars:**

Coloured bars are included for reference when instructing the patient - For example, "Start with the line below the red bars". The green bars also indicate the normal driving license standard.



#### **Adults Letter Charts**

◆ Sizes: 6/60, 6/36, 6/30, 6/24, 6/18, 6/12, 6/9, 6/6, (6/5)

◆ Printer: Colour printer - laser or ink-jet

◆ Paper: ● Standard A4 paper
For use at 3 metres (10 feet)

◆ Paper: **②** A3 paper - 150% enlargement For use at 6 metres (20 feet)

◆ Note: This is the most popular chart. It is provided in three versions to counter memorising:

◆ Ref E - Top letter E.

◆ Ref D - Top letter D.

◆ Ref T - Top letter T.

◆ Ref E5 - Top letter E - Incl 6/5.



#### Adults Letter Charts - Indirect

◆ Sizes: 6/60, 6/36, 6/30, 6/24, 6/18, 6/12, 6/9, 6/6, (6/5)

◆ Printer: Colour printer - laser or ink-jet

◆ Paper: ● Standard A4 paper For use at 3 metres (10 feet)

◆ Paper: **②** A3 paper - 150% enlargement For use at 6 metres (20 feet)


Note: This chart is designed for viewing through a mirror ("indirectly"). It is provided in two versions to counter memorising:


◆ Ref E - Top letter E.


◆ Ref D - Top letter D.

◆ Ref E5 - Top letter E - Incl 6/5.

Document: Info\_V\_01\_Descrip Last edited: 29/12/2022 Author: Anthony Hanks, OD







#### **Adults Number Charts**

◆ Sizes: 6/60, 6/36, 6/30, 6/24, 6/18, 6/12, 6/9, 6/6, (6/5)

◆ Printer: Colour printer - laser or ink-jet

◆ Paper: **①** Standard A4 paper For use at 3 metres (10 feet)

◆ Paper: **②** A3 paper - 150% enlargement For use at 6 metres (20 feet)

♦ Note: These charts are useful for patients who do not recognise letter, or as an alternative chart. It is provided in two versions to counter memorising:

◆ Ref N3 - Top number 3.

◆ Ref N8 - Top number 8.

◆ Ref N35 - Top number 3 - Incl 6/5.

## **Adults LogMAR Charts**

◆ Sizes: 0.8, 0.7, 0.6, 0.5, 0.4, 0.3, 0.2, 0.1, 0.0, -0.1, -0.2

◆ Printer: Colour printer - laser or ink-jet

◆ Paper: ● Standard A4 paper For use at 3 metres (10 feet)

◆ Paper: **②** A3 paper - 150% enlargement For use at 6 metres (20 feet)

◆ Note: These charts are based upon the Logarithm of the Minimum Angle of Resolution (LogMAR) and are often used in vision research. It is provided in two versions to counter memorising:

Ref C - Top letter C.

Ref H - Top letter H.

 $LogMAR = Log_{10} (1 / Decimal VA)$ 

## **Children's Optotype Charts**

◆ Sizes: 6/60, 6/36, 6/30, 6/24, 6/18, 6/12, 6/9, 6/6, (6/5)

◆ Printer: Colour printer - laser or ink-jet

◆ Paper: **①** Standard A4 paper For use at 3 metres (10 feet)

◆ Paper: **②** A3 paper - 150% enlargement For use at 6 metres (20 feet)

◆ Note: For use with young children who cannot recognise letters reliably.


The symbols (*Hanks Optotypes*) have been tested for recognition with preschool children (aged 3 to 5 years).


This chart is provided in two versions to counter memorising:


◆ Ref H - Top symbol Heart.

◆ Ref A - Top symbol Apple.

◆ Ref H5 - Top symbol Heart - Incl 6/5.







## **Tumbling E Charts (Illiterate E)**

Sizes: 6/60, 6/36, 6/30, 6/24, 6/18, 6/12, 6/9, 6/6, (6/5)

Printer: Colour printer - laser or ink-jet

Paper: • Standard A4 paper For use at 3 metres (10 feet)

Paper: 2 A3 paper - 150% enlargement For use at 6 metres (20 feet)

Note: For use with young children, people who cannot speak English and others who cannot recognise letters.

> This chart is provided in two versions to counter memorising:

◆ Ref R - Top symbol Right.

◆ Ref D - Top symbol Down.

◆ Ref R5 - Top symbol Right - Incl 6/5.

### **Landolt Ring Chart**

Sizes: 6/60, 6/36, 6/30, 6/24, 6/18, 6/12, 6/9, 6/6, (6/5)

Printer: Colour printer - laser or ink-jet

Paper: • Standard A4 paper For use at 3 metres (10 feet)

Paper: 2 A3 paper - 150% enlargement For use at 6 metres (20 feet)

Note: For use with young children, people who cannot speak English and others who cannot recognise letters.

This chart is provided in one version:

◆ Ref C - Top symbol C.

◆ Ref C - Top symbol C5 - Incl 6/5.

## **Astigmatism Charts**

Printer: Colour printer - laser or ink-jet

Paper: A4 or A3 paper

3 to 6 metres (10 to 20 feet)

Note: This chart is used when screening

for astigmatism, or measuring the cylinder axis of the astigmatism.

There are two astigmatism charts: ◆ Astigmatism Clock & Blocks.

Astigmatism Axis Chart.

## Notes About Astigmatism Interpretation

The red numbers on the Astigmatism Clock Chart are desgigned for the patient, rather than for the clinician. They do NOT indicate the axis directly - for example 11 on the chart does not indicate an axis of 110.

This is descibed in our *Eyecare Education Series* E-03 on page 32 of this document.